

IE Carefully consult the instructions for use before using Mycosan Athlete's Foot Precision brush, even if you have used this product before. The information may have been amended according to new research and findings since you last used the product. Store these instructions for use for future reference.



Mycosan Athlete's Foot Precision brush

Suitable for diabetics and children of 4 years and older.



SUITABLE FOR
DIABETICS

Mycosan Athlete's Foot Precision Brush contains 15 ml gel in a precision brush applicator tube.

Athlete's foot is a fungal infection of the skin of the foot.

Mycosan Athlete's Foot is a hygienic and effective precision treatment. It is indicated for treatment and relief of symptoms of athlete's foot by providing an unfavourable environment for fungal growth by acidification.

Treatment is important to prevent the fungal infection from spreading or worsening.

Why use Mycosan Athlete's Foot?

- ✓ **Precision Brush:** hygienic & effective.
- ✓ **Clears athlete's foot symptoms:** e.g. itching, burning, cracking, redness, flaking of skin and bad odour.
- ✓ **Restores:** skin is regenerated.
- ✓ **Suitable:** for diabetics and children 4 years+.

How does the Mycosan precision brush work?



1. HYGIENIC:
For fast and hygienic application.



2. REACHES FUNGI BETWEEN TOES:
For targeted application between the toes to optimize result.

Mycosan Athlete's Foot has a unique two-fold effect:

1. Instant itch relief.
2. Restores the skin.

Dosage and use

Apply a thin layer twice a day to the affected area(s) as early as possible following the onset of signs and/or symptoms, i.e. itching, cracking, burning, redness, malodour and flaking of the skin. Treatment should be continued until healing is complete. Continue treatment for 2 weeks after the symptoms have cleared to prevent recurrence. Consult a health care professional if the symptoms have not improved after 2 weeks.

Ingredients

Triacetin (40%), PEG-6, alcohol denat., glyceryl undecylenate, aqua, dimethyl Isosorbide, rye ferment filtrate, pentylene glycol, hydroxypropylcellulose, parfum, cetylpyridinium chloride, terpineol, terpinolene, turpentine, limonene, alpha-terpinene, pinene, citral, eucalyptus globulus oil, camphor.

What is athlete's foot?

Athlete's foot or tinea pedis is a fungal infection of the skin of the feet. Three different forms can be distinguished:

Interdigital tinea pedis is an infection of the web spaces, most common between the fourth and fifth toe, and spreading to the third or further. It itches, the skin appears white and macerated, with weepy red erosions beneath, and it is often malodorous because of a bacterial superinfection.

Moccasin-type tinea pedis presents as bran-like scale, with silvery white flakes on a red thickened base distributed over the sole, heel, and sides of the feet, and tends to be chronic.

Vesicular tinea pedis usually occurs on the sole, but also on the side of the foot. Vesicles or pustules are present and the area may become macerated and secondarily infected with bacteria, which causes an unpleasant odour.

Without treatment, the condition is likely to become worse and the infection may spread to other areas, the nails, or even contaminate other individuals.

Many factors promote athlete's foot: warmth, occlusive shoes, certain sports including swimming, hiking or marathon running and, generally speaking, all situations in which feet are prone to sweating, maceration and dampness.

Athlete's foot may also result from fungal contamination after prior infection of the nail, as for example in onychomycosis (fungal nail). In fact, when poorly treated or left entirely untreated, nail fungal infections can spread and cause athlete's foot. It should be noted that certain individuals are more easily subject to fungal foot infections: diabetics, people suffering from blood circulation problems, and those with occupations requiring the wearing of safety shoes or boots, etc.

Use in combination with other products

When using a foot cream, apply Mycosan Athlete's Foot at least 15 minutes in advance.

Precautions for use

- For personal use only and should not be shared with others.
- Suitable for local external use on the skin only.
- Do not use on open wounds.
- Do not use if allergic or oversensitive to one or more of the ingredients. Stop the use immediately when irritation or an allergic reaction occurs. Contact a health care professional if symptoms persist.
- When pregnant or breastfeeding, only use if instructed by and under supervision of a health care professional.
- Keep out of reach of children.
- Squeeze tube gently and close the tube securely after use.
- It is recommended to consult a physician when treating large areas for an extended period.

Recommendations to optimise the efficacy of Mycosan Athlete's Foot:

- Dry your feet after washing, including between the toes.
- Opt for wearing absorbent socks to keep your feet fresh and dry. Change them regularly in case of high perspiration.
- Change shoes regularly and let them air and dry out if damp. Avoid wearing shoes made from synthetic materials.
- Use an absorbent protective powder on feet and in shoes to keep the feet dry.
- Choose well-fitting shoes that are not too narrow or too tight; allowing your feet to breathe properly.
- Always wear footwear in high-risk public places (changing rooms, pools, saunas etc.).

Possible side effects

Mycosan Athlete's Foot may sting a little when applied to breached skin.

If a serious incident has occurred in relation to the product, report to the manufacturer mentioned on the leaflet and the competent authority in your country.

FOR MORE INFORMATION, PLEASE VISIT WWW.MYCOSAN.COM



Shelf life and storage

Store at room temperature.
Keep out of reach of children.



Keep tube away from sunlight.



Single patient, multiple use



Lot number and expiration date: see box & tube

When the treatment is complete dispose of the tube, leaflet and carton as with any common urban waste. Consider the environment and recycle.



Manufacturer:

Serrix B.V.
Herengracht 458
1017 CA Amsterdam
The Netherlands



Distributor:

Pamex Limited
Moneen Business
Park, Castlebar, Co.
Mayo, F23 RX24,
Ireland
094 9024000
ask@pamex.ie



Medical Device